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Modelle

Typ G90 – SUPREME line:

**79-130660-xyz(L), 79-200660-xyz(L),
79-200960-xyz, 79-201260-xyz, 79-201460-xyz**

Typ G30 – ECO line:

**73-200660-xyz(L), 73-200960-xyz,
73-201260-xyz, 73-201460-xyz**

*Operating Instructions
Instructions de service
Bedieningshandleiding
Istruzioni per l'uso
Instrucciones de servicio*

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Type G30 – ECO line:

73-200660-xyz(L), 73-200960-xyz, 73-201260-xyz, 73-201460-xyz

Please enter the following details:

(this information is important in the event of enquiries)

Cabinet model: _____ Production number: _____

Key number: _____ Location: _____

Responsible person: _____

Operating instructions CONTENTS

Section		Page
1	Manufacturer and technical details of the safety cabinet and safety instructions	15
2	Transport of the safety cabinet	16
3	Installing the safety cabinet	17
4	Putting the safety cabinet into service	17
5	Interior fittings of the safety cabinet	17
6	Venting	17
7	Opening the safety cabinet after a fire	19

Section		Page
8	Testing/fault	19
9	Disposal	20
10	Guarantee	20
11	Spare parts	20
12	Assembly instructions	9
13	EC Declaration of Conformity	55
14	GS Certificate	56
15	All-round service	21

1. Manufacturer and technical details

1.1 Sales and customer service

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Mainparkstrasse 6-10
63801 Kleinostheim
Germany

Phone +49 6027 403-0
Fax +49 6027 403-121

E-mail: info@dueperthal.com
www.dueperthal.com

1.2 Name/application

Safety cabinet according to DIN EN 14470-2 – Type G90/Type G30 with a maximum fire resistance of 90 minutes/30 minutes respectively and a temperature increase of 50 Kelvin measured at the neck of the gas cylinder.

1.3 Construction

Outer carcass: Powder coated sheet steel.
Insulation: Insulation panels with sandwich construction.
Inner carcass: Melamine resin coated decor panels, cabinet inside floor made of sheet steel.

Code for x in the Ref. No.:

0 = Carcass light grey RAL 7035, door(s) ruby red RAL 3003.
4 = Carcass light grey RAL 7035, door(s) light grey RAL 7035.
8 = Carcass and door(s) in special colour.
9 = Carcass light grey RAL 7035, door(s) in special colour.

The production number (serial number) is printed on the adhesive label, e.g. 40001.

1.4 Fittings and features

Model type	Model size				
	S	M	L	XL	XXL
Type G90 – SUPREME line	79-130660-xyz	79-200660-xyz	79-200960-xyz	79-201260-xyz	79-201460-xyz
Type G30 - ECO line	-	73-200660-xyz	73-200960-xyz	73-201260-xyz	73-201460-xyz

(Variable code: x = colour, yz = definition of the interior fittings/features, L = optional door hung left, R = interior fittings, reversed)

1.5 General technical details

Model	Type G90 – SUPREME line					Type G30 – ECO line			
	S	M	L	XL	XXL	M	L	XL	XXL
Outer dimensions:									
Width approx. mm	595	595	895	1195	1395	595	895	1195	1395
Depth approx. mm	610	610	610	610	610	595	595	595	595
Height approx. mm	1444	2099	2099	2099	2099	2099	2099	2099	2099
Empty weight approx. kg	290	390	530	660	740	310	440	523	655
Max. total weight when fully loaded approx. kg	390	490	710	920	1080	410	620	783	995
Volume in approx. m ³	0,22	0,33	0,56	0,78	0,95	0,40	0,64	0,88	0,98
Standing spaces for 50 l cylinders	–	1	2	3	4	1	2	3	4
Load bearing capacity per intermediate shelf (optional) with uniformly distributed load in approx. kg	50	50	50	50	50	50	50	50	50

1.6 Guarantees

Our liability for delivery defects is specified in our terms and conditions of delivery. We shall not accept any liability whatsoever for damage caused by failure to observe the operating and maintenance instructions and the conditions of use.

Only faultless products which have the properties assured by us leave our factory. Thus, if the following information provided in these operating and maintenance instructions is observed, the product should work properly.

1.7 Safety

These operating and maintenance instructions contain fundamental information to be observed for the operation and maintenance/repair. Therefore, these operating and maintenance instructions shall always read by the owner/operator's skilled personnel before putting the product into service and the owner/operator shall draw up their own in-house instructions for use based on these operating instructions and the local and specific operational circumstances and shall instruct the responsible skilled personnel accordingly before they start working with the product.

The safety information contained in these operating instructions which, if not observed, can cause risks to people and the environment, is specially labelled with general warning symbols in accordance with the accident prevention regulations BGV A 8

to warn of electrical voltage with.

In case of safety instructions which, if ignored can cause risks to the technical equipment and its functions,

the word **IMPORTANT** is added.

Notices directly attached to the safety cabinet must always be observed and maintained in a fully legible condition.

1.8 Risks in the event of failure to observe the safety instructions

Failure to observe the safety instructions can result both in risk to people as well as the environment and technical equipment. Failure to observe the safety instructions can result in all guarantees and compensation claims becoming null and void.

The safety instructions given in these operating and maintenance instructions, the existing legal regulations as well as the accident prevention regulations and the owner/operator's company and safety regulations shall be observed.

Safety instructions for the owner/operator

- the bans on storing hazardous substances and preparations together are to be observed,

- the owner/operator shall forbid unauthorised persons to access the hazardous substances stored in the safety cabinet,
- the owner/operator shall regulate smoking, handling and working with naked flames and spark-drawing work in and on the safety cabinet,
- safety devices are to be kept fully functional at all times,
- the function of the inlet and air extraction openings of the industrial venting shall not be impaired,

Safety instructions for maintenance, inspection and installation work

The owner/operator shall ensure that all maintenance and installation work is carried out by authorised, qualified personnel.

All work on the electrical installations shall only be carried out by qualified electricians when the system is off-load – see also the relevant accident prevention regulations, the VDE regulations and the provisions of the local power supply company.

Impermissible use

The safety devices and features installed in the DÜPERTHAL safety cabinet, e.g. industrial venting, shall be in full working order and kept so at all times.

Impermissible use is also deemed to be granting unqualified or inadequately personnel and unauthorised persons free access to the hazardous substances stored in the safety cabinet.

1.9 Possible uses

DÜPERTHAL safety cabinets are approved for the storage of pressurised gas cylinders in working spaces. The internal height of the safety cabinet ensures that standard commercial pressurised gas cylinders including the necessary pressurised gas fittings can be placed in the cabinet and fitted. By inserting shelves for side installation, smaller cylinders, e.g. 5 and 10 litre, can be flexibly accommodated.

NOTE

The volume of all pressurised gas cylinders in safety cabinets may not exceed 220 litres, this includes pressurised gas cylinders for purging gases with a maximum volume of 10 litres. Depending on the model, up to 4 pressurised gas cylinders size 50 litre or up to 3 pressurised gas cylinders size 70 litre can be stored in a safety cabinet.

2. Transport

2.1 Transport safety measures

The safety cabinet is equipped with a transport skid to protect the cabinet during transport. Further, an industrial truck, e.g. hand pallet truck 1.31.005 from DÜPERTHAL can run under the safety cabinet. Please also note item 2.2.

2.2 Safety measures for in-house transport

The unloaded safety cabinet may only be picked up from the underside using suitable means of transport and must be

transported upright. It is very important to avoid damage due to tilting when lifting the cabinet for transport. Tilting can cause damage, which restricts or destroys the functions of the installed safety technology.

Please note and follow the assembly instructions on page 9.

3. Installation

3.1 Installation location requirements

The surface must be level and flat and be able to bear the weight of the safety cabinet when fully loaded.

Furthermore, the installation location must ensure that the safety cabinet:

- cannot be damaged by vehicles.
- is not exposed to any direct or indirect heat source.
- is protected against moisture.
- is not exposed to an operating temperature below -5 °C or above +45 °C.

3.2. Installation

Safety cabinets must be installed horizontal. Only by aligning the safety cabinet do you ensure that it cannot overturn and also, perfect function of the safety features is only ensured under this condition.

4. Putting into service

4.1 Wing door(s)

The lever handle must be turned through at least 90° to open the wing door. The wing door can then be pulled open using the lever handle. To open the left-hand wing door of double-wing doors, the locking device on the end face must also be rotated through 90°. This unlocks the wing door and enables it to be pulled open. The doors must be manually closed in the reverse order after each use.

4.1.1 Lock cylinder

The door can be locked using the integrated lock cylinder. The lock cylinders can be subsequently adapted to the user's needs.

4.2 Storing different types of gas together

Under state of the art standards (technical rules for pressurised gases, TRG 280) it is permitted to store different gases in a safety cabinet. Gases which can react with each other must be stored separately. Please also read and follow the EC safety data sheets for the relevant gases.

4.3 Storing corrosive gases

The storage of corrosive gases can have effects on the function of the supply and exhaust air shut-off devices. DÜPERTHAL does not accept any responsibility whatsoever for damage caused by corrosion.

4.4 Earthing possibility (optional)

An earthing device inside the cabinet is possible as an optional feature. The equipotential bonding saddle is already provided on the cabinet's roof.

5. Interior fittings

5.1 Pressurised gas fittings

Installation rails are provided on the inside rear panel of the cabinet for fixing pressurised gas fittings. The fittings can be screwed directly onto the installation rails.

5.2. Installation and use of cylinder holders

The continuous cylinder holder is fixed to the side panels in the factory. Individual height adjustment is possible. Undo the screws, align the cylinder holder at the required height and screw onto the side

panels. Each cylinder standing space is equipped with a tension belt for fixing the pressurised gas cylinders.

5.3 Shelf for side installation

Attach the side shelf at the required height on the right or left using the screws provided. The cylinder holder including tension belt is included in the scope of supply. The load-bearing capacity is 50 kg for uniformly distributed load.

5.4 Pipe penetration (for pressurised gas cylinders in use)

The number of pipe penetrations must be limited to the minimum number required and may not exceed three per cylinder. The predefined positions on the roof of the safety cabinet are to be used for the pipe penetrations. The side and rear panels are not approved for pipe penetrations. The maximum allowable diameter per pipe is 10 mm and in general may not exceed the nominal diameter of the fed through pipe.

Penetrations and holes no longer used must be properly closed off – please contact DÜPERTHAL.

IMPORTANT It is advisable to make the pipes from stainless steel or a material with similar thermal conductivity, to maintain the fire resistance according to which the safety cabinets are classified.

5.5 Installing electric cables

The number of cable penetrations must be limited to the minimum number required and may not exceed two per cylinder. The predefined positions on the roof of the safety cabinet are to be used for the cable penetrations. The side and rear panels are not approved for pipe penetrations. The maximum allowable diameter per cable is 20 mm and in general may not exceed the nominal diameter of the fed through cable.

Penetrations and holes no longer used must be properly closed off – please contact DÜPERTHAL.

5.6 Roll-in flap

The roll-in flap can be easily folded down by foot operation. When pressed the flap slowly lowers, pneumatically damped (without external energy). After placing the pressurised gas cylinder in the cabinet, and if necessary connecting it, the flap must be folded up again and the doors closed.

6. Venting

6.1 General information

According to DIN EN 14470-2, the following air exchange is specified for a safety cabinet connected to an industrial exhaust air system, with the cabinet doors closed and locked:

- at least 10-fold air exchange/hour for flammable and oxidising gases.
- harmful, toxic and highly toxic gases at least 120-fold air exchange/hour.

The pressure drop in the safety cabinet may not exceed 150 Pa under the above prerequisites. The venting system must cause a partial vacuum in the cabinet. Both the top and bottom half of the cabinet must be vented. The exhaust stream must ensure that small quantities of leaking gases are extracted. Extraction can take place via a separate ventilator or via a special extractor and must be in operation 24 hours a day (VDI 2051, Item 3). The exhaust stream must ensure that small quantities of leaking gases are extracted. The exhaust air system must lead to a safe place in outdoors, i.e. discharge to the open air. In the event of a fire, self-closing supply and exhaust air valves prevent heat from penetrating.

Checking of the venting must include a visual inspection as well as measurement of the air stream and the pressure drop inside the empty cabinet.

NOTE The storage of corrosive gases can have effects on the function of the supply and exhaust air opening shut-off devices.

6.1.1 Pressure drop in the safety cabinet

SUPREME line - Type G90				
Model	10-fold air exchange/h		120-fold air exchange/h	
	Volumetric flow rate in m ³ /h	Pressure drop in the cabinet	Volumetric flow rate in m ³ /h	Pressure drop in the cabinet
79-130660-xyz	2,2	< 1 Pa.	27	< 150 Pa.
79-200660-xyz	3.3	< 1 Pa.	40	< 150 Pa.
79-200960-xyz	5.6	< 1 Pa.	67	< 150 Pa.
79-201260-xyz	7,8	< 1 Pa.	94	< 150 Pa.
79-201460-xyz	9.5	< 1 Pa.	114	< 150 Pa.

ECO line - Type G30				
Model	10-fold air exchange/h		120-fold air exchange/h	
	Volumetric flow rate in m ³ /h	Pressure drop in the cabinet	Volumetric flow rate in m ³ /h	Pressure drop in the cabinet
73-200660-xyz	4.0	< 1 Pa.	54	< 150 Pa.
73-200960-xyz	6.4	< 1 Pa.	77	< 150 Pa.
73-201260-xyz	8.8	< 1 Pa.	106	< 150 Pa.
73-201460-xyz	9,8	< 1 Pa.	118	< 150 Pa.




6.2 Connection to an exhaust air system

The exhaust air and air supply openings are located on the top of the cabinet of the safety cabinet. The scope of supply includes one exhaust air and one air supply connection NW Ø 75, which enables connection to a ventilation system, e.g. DÜPERTHAL 2.00.320 for exhaust air monitoring with ventilator (ATEX-compliant).


The venting pipe must be connected to the exhaust air connection socket - left-hand socket viewed from the front. The pipe must be connected to the exhaust air connection socket using a collar, or similar. The air supply can be taken from the room and under normal operating conditions does not require any additional air supply from the outside. After installing the safety cabinet it must be checked to ensure it has been properly connected to the extractor; this includes a visual inspection as well as measurement of the air stream and the pressure drop inside the empty cabinet. The installation of an industrial venting system or connection to an existing exhaust air system is not part of the DÜPERTHAL scope of supply. Please contact a local ventilation company near you.

IMPORTANT The attachment of a ventilator represents an equipment combination. As Zone 2 equipment, the ventilator must fulfil the requirements of the EC Directive 94/9/EC (ATEX) equipment group II, category 3 G.

6.2.1 Top views

(Variable code: x = colour, yz = definition of the interior fittings/features, L = optional door hung left, R = interior fittings, reversed)

6.3 Air ducting system inside the cabinet

The slits in the air ducts make the venting effective in both the top and bottom half of the safety cabinet.

6.4 Choice of ventilator

Only ventilators which fulfil the technical specifications defined in the VDMA leaflet 24 169 Part 1 may be used. In non-hazardous areas this means according to VDMA, Indoors: Zone 2 - outdoors: non-hazardous area.

In hazardous areas (i.e. potentially explosive atmospheres) the motor must also be explosion-proof.

IMPORTANT

The attachment of a ventilator represents an equipment combination. As Zone 2 equipment, the ventilator must fulfil the requirements of the EC Directive 94/9/EC (ATEX) equipment group II, category 3 G. See also the DUPERTHAL catalogue PA 4 "Ventilators".

7. Opening the safety cabinet after a fire

After a fire the safety cabinet may not be opened for at least 24 hours after the fire.

IMPORTANT

Depending on the fire duration and load, an explosive gas - air mixture may have formed inside the safety cabinet. Before opening the safety cabinet, all ignition sources, naked light, etc. within a perimeter of 10 m around the safety cabinet shall be removed. If the safety cabinet has to be opened with tools, spark free tools only must be used.

8. Testing/Fault/Markings

The safety cabinet is to be checked for externally identifiable damage or defects after installation and before putting into service and thereafter at the following intervals as well as after changes or repairs.

8.1 Monthly check

The monthly check includes the following items:

- Door closure: Open the door and check its closing.
- Venting: Test the effectiveness of the venting using a woollen thread, smoke tube or similar at the slits of the exhaust air venting duct.
- Sealing strips: Check the sealing strips correctly fit in the door rebate and the end face(s) of the door(s). If the sealing strips are visibly damaged they must be immediately replaced. Please contact the service/customer service department.

8.2 Annual check

The annual check (inspection) of the safety cabinet must be carried out by a competent, qualified person. If faults occur you can help the technical customer service (see cover for address and phone number) by quoting the cabinet model and production number and by describing the fault.

8.3 Markings

The following must be attached in a suitable place on the front of the safety cabinet so that it is clearly visible:

- Fire resistance in minutes, e.g. Type G90
- Notice: "Close door"
- Pictogram or sign according to ISO 3864 with respect to the pressurised gas cylinders stored
- Name and/or trademark of the manufacturer
- Serial number and year of construction
- Label the exhaust air and air supply connection sockets to avoid mix-ups
- Note instructing user to read the operating instructions.

The size of all symbols and information signs should be appropriate for the size of the safety cabinet.

Please clean cabinet stickers with normal commercial glass cleaners only. Other cleaners may attack the sticker.

Possible additional markings/labels (installed by customer on site):

- Warning symbol W01 "Warning of flammable substances" and/or warning symbol W03 "Warning of toxic substances".
- Note that the safety cabinet is only intended for the storage of inert gases.

8.4 Faults

Faults	Possible cause	Possible solutions
Door cannot be closed.	Safety cabinet is not plumb.	Install the safety cabinet horizontally (align).
No air exchange.	Shut-off device is closed, as fuse cutout is defective.	Fuse cutout to be replaced by technical customer service/service.

9. Disposal

The safety cabinet can be completely dismantled. The individual fractions, e.g. metal, insulation panels, etc. can be separately recycled or reused. The national and local disposal regulations shall be observed. To conserve resources, parts of the safety cabinet or the whole safety cabinet should not be thrown away in the bulk or household waste.

10. Guarantee

All damage caused by improper use and installation and corrosive damage are excluded from the guarantee.

11. Spare parts/accessories

- Product:**
- Sealing strips
 - Shut-off device
 - Handle
 - Installation rails
 - Intermediate shelf
 - Cylinder holder
 - Lock cylinder

General information

These operating instructions take into account the requirements of the relevant legal provisions as of May 2004.

IMPORTANT

These operating instructions will not be updated by us in the event of changes to the relevant laws, regulations, etc., unless we specifically asked to update them.
These operating instructions do not replace the company work instructions and operating instructions to be prepared by the owner/operator and required by law **for safety cabinet for pressurised gas cylinders.**

All dimension are approximate.
Subject to technical amendments. Printing errors excepted.

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All rights reserved. Reprint prohibited, even as extracts. No part of this document may be reproduced or processed using electronic systems, copied or distributed in any form whatsoever (photocopy, microfilm or other method, not even for teaching purposes, without the prior consent of DÜPERTHAL.

Thank you for purchasing

DÜPERTHAL SAFETY CABINETS.

To ensure the safety function of your safety cabinet you should have it checked by our experienced and competent service technicians once a year.

Our tip: **All-round service**

with the DÜPERTHAL Service Agreement **DS 2/5-4***

2
You secure
2 decisive
advantages:

5
With a term of 5 years
the manufacturer's guarantee** is extended to
5 years.

-4
The costs for the first
inspection are borne by
us, you only pay for
4 years.

All-round service at a glance:

- Up to 5 year manufacturer's guarantee.**
- First inspection free.
- Cost savings for repairs, time is not billed separately (included in lump sum charge).
- **Travel costs are included in the lump sum charges.**
- Your cabinet's protective function is retained permanently.
- Support from trained and certified service technicians.
- No tiresome scheduling, we automatically remind you of pending inspection dates.
- **DS 2/5-4*** – Customers are given preferential treatment in the scheduling.

All-round service

Customer service

Phone
+49 6027 403-157
service@dueperthal.com

www.dueperthal.com

DÜPERTHAL is certified to
DIN EN ISO 9001

Client:

(please state precise address and phone no.)

Service agreement No.:

between the recipient as
 Client and
**DÜPERTHAL SICHERHEITSTECHNIK
 GmbH & Co. KG**
 Mainparkstrasse 6-10
 63801 Kleinostheim
 Germany

The annual inspection of safety cabinets with proof of testing is legally required by the Arbeitsstättenverordnung (workplace regulations). Employers are responsible for implementation.

The DÜPERTHAL Service Agreement DS 2/5-4* includes the following services:

- Annual safety inspection and maintenance
- Test certificate and test sticker
- Time required for repairs (spare parts are billed separately)
- **Travel costs**
- 5 year manufacturer's guarantee**

Number of cabinets:

1 _____ (please enter quantity)
 2 - 3 _____
 4 - 6 _____
 7 - 10 _____
 11 - 20 _____
 20 plus _____

Lump sum per cabinet/year:

The inspection and maintenance covers all **safety relevant parts**, such as fire protection and air seals, fusible links, ventilation openings, door closers and locks and the insulation. Furthermore, all **mechanical components** are checked to ensure that they are in sound working order and an **exhaust air measurement is taken**. Malfunctions are immediately corrected.

Deviations from this contract are only valid if they have been agreed in writing by both parties to the contract. Otherwise the Service Terms and Conditions (supplement) shall be complied with.

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 (Stamp/Signature)

 (Stamp/Signature)

* Only valid for cabinets less than 2 years old (from delivery date).
 ** Manufacturer's guarantee does not apply to third party makes.

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EG-DECLARATION OF CONFORMITY

im accordance with Machinery Directive 2006/42/EG, Appendix II A

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hereby declare, that the design, construction and commercialized execution of the below mentioned machine complies with the essential health and safety requirements of the EG Machinery Directive. The validity of this declaration will cease in case of any modification or supplement not being agreed with us previously. Furthermore, validity of this declaration will cease in case that the machine will not be operated correctly and in accordance to the operating instructions and/or not to be inspected regularly.

Machine description:

safety cabinet

Machine type:

73-201260-xyz / 73-200960-xyz / 73-200660-xyz(L) /
73-200670-xyz(L) / 73-130660-xyz(L)79-201260-xyz / 79-200960-xyz / 79-200660-xyz(L) /
79-200670-xyz(L) / 79-130660-xyz(L)7w-201260-xyz / 7w-200960-xyz / 7w-200660-xyz(L) /
7w-200670-xyz(L) / 7w-130660-xyz(L)variables:

w = Type G - class

x = colour

yz = definition of the interior

L = optional door hang left

Relevant EG Directives:

Machinery Directive (2006/42/EG)Transposed national standards
and technical specifications in
particular:

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Parameters:

key code wxyz

w stands for number 1, 3, 6 or 9, description type of
fire resistance according DIN EN 14470-2:2006

1 = Typ G 15, 3 = Typ G 30

6 = Typ G 60, 9 = Typ G 90

x colour description, stands for number 0, 4, 5, 8, 9

yz description of interior fitting,
stands for number 00 to 99

In case of door hinge left side

L will be added at the article-number.

Tested according to: DIN EN 14470-2:2006
DIN EN 14727:2006

The product meets the requirements of the German Equipment and Product Safety Act. The Certification marks shown above can be affixed on the product. The certification marks must not be altered in any way. The use of the GS-Mark is permitted until the listed date, the use of the TÜV-Mark is unlimited, unless it is cancelled. See also notes overleaf.

Test report no.: 71323438

GS-Mark valid until: 2012-06-22

Date, 2007-06-22

Page 1 of 1

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© Copyright 2006